





Introducing XI-S+™

For Soft Tissue Repair



+ BRINGING A
NEW PERSPECTIVE
TO PATIENT CARE



www.co-therapeutics.com



Introducing the Colorado Therapeutics XI-S+™ For Soft Tissue Repair

Delivering a biocompatibility and
durability profile to address patient needs.

XI-S+™ is a xenogenic (porcine) pericardium that undergoes a proprietary Colorado Therapeutics cross-linking tissue treatment process that provides a strong, durable, and biocompatible product with extremely low glutaraldehyde and DNA residuals. The unique benefits of XI-S+™ are well-suited for repairing hernia defects and for use in applications where current xenograft soft tissue repair products are indicated.

THE COLORADO THERAPEUTICS XI-S+™ has received FDA 510(k) clearance and is intended for use as a soft tissue implant to reinforce soft tissue where weakness exists and for the surgical repair of damaged or ruptured soft tissue membranes.

Colorado Therapeutics XI-S+™ is for single use only and cannot be resterilized.

XI-S+™ is the first in a family of xenologic products from Colorado Therapeutics LLC. Colorado Therapeutics has a deep research pipeline of products under development to address unmet patient and market needs in soft tissue repair.

● ORDERING INFORMATION:

PRODUCT NUMBER	PRODUCT SIZE	PRODUCT NUMBER	PRODUCT SIZE
X10203-S+	2 cm x 3 cm	XI0710-S+	7 cm x 10 cm
X10404-S+	4 cm x 4 cm	XI0808-S+	8 cm x 8 cm
XI0606-S+	6 cm x 6 cm	XI0615-S+	6 cm x 15 cm
XI0508-S+	5 cm x 8 cm		

Contact your sales representative or Colorado Therapeutics LLC for ordering information, in-service support materials, and technical specifications.

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician. For complete instructions for use, storage, warnings, indications, contraindications, precautions, adverse reactions, and disclaimer of warranties, please refer to the insert accompanying each Colorado Therapeutics Xenograft Implant product.

Colorado Therapeutics products are protected by trademarks and know-how. 01/18 © Colorado Therapeutics LLC. PK038.02 CT. All rights reserved. Printed in the U.S.A.

Colorado Therapeutics LLC
2150 West 6th Avenue, Ste L
Broomfield, CO 80020

877.290.0045 (toll-free)
303.469.9459
720.699.1391 (fax)

orders@co-therapeutics.com
www.co-therapeutics.com

